

Contents

1	Power-On Procedures (IMPORTANT!)	2
2	Project Organization	2
3	Data Format	3
3.1	Saving and Loading	3
3.2	Editing Files	4
4	Settings	4
5	Characterizing MZMs	4
6	Characterizing Lasers	4
7	Characterizing Symetrix 402 Delay Module	5
8	Characterizing bandpass filter	6
9	The Circuit	6
9.1	The AD5206 Digital Potentiometer	6
10	Utilities	10
10.1	combine_figs.py	10

1 Power-On Procedures (IMPORTANT!)

To avoid damaging the DigiPot chip or the Arduino, circuit components must be turned on in a very particular order:

1. Turn on the power supply to the DigiPot ($\pm 2.5\text{V}$). If this is the first time the circuit has been turned on in a while, check that the voltages are reasonable with a multimeter.
2. Turn on Arduino power by plugging in the USB cable.
3. Power on delay module to ensure grounding.
4. Turn on opamp power ($\pm 15\text{V}$) using the breadboard switch.
5. Power lasers up (typically, to about 20mA).
6. Turn on photodetectors—by powering on the lasers before the detectors, we can minimize the effect of the relaxation oscillation and associated current surge.
7. Ready

2 Project Organization

Figure 1 shows how the project is organized.

- The data directory contains the raw data collected.
 - delay contains data on the behavior and performance of the delay modules
 - digipot contains information on the AD5206 digital potentiometer (a chip you probably want to avoid)
 - filters has information on the bandpass filters and their characteristics over various frequency ranges
 - lasers has current-power curves for the three lasers in Lucas' lab.
 - machzehnder has information on the response of the three MZMs.
- doc has this document, and associated source files.
- figures has various figures and results, derived from the data. It has the same categories as the data directory.
 - The pinouts directory has pinouts for the various ICs I used.
- manuals has PDF manuals for the various bits used in the circuit.
- notes has application notes and other writeups. These also have copies in the lab notebook.
- numericalanalysis has my work on numerical bifurcation and theoretical prediction of the size and shape of the amplitude death region.
- papers has PDF copies of various relevant papers.
- scripts has all the Python source code for collecting and analyzing data. I'll get into more depth on what each of these scripts do later.




Figure 1: The directory format of the project.

3 Data Format

All data is stored in the Numpy builtin compressed-archive format (.npz), since it's easy to work with in Python, makes it possible to store the data and metadata about its capture (such as time, test parameters, etc.) in the same file, and is nice and compact (no more 80MB captures!).

3.1 Saving and Loading

I've written a few utility functions to make saving and loading these files as simple as possible. Located in `utils.py`, the `save` function takes a filename, arbitrary Numpy array (any shape or dimension number), and (optionally) a dictionary of settings to save; and writes them to a file.

Similarly, the `load` function, given a filename, loads that file and returns a `data`, `metadata` tuple. If the metadata is stored in a separate file (which I was doing for the first few weeks), this function loads from that instead, and resaves the data in compressed form. The code snippet below shows how to use these functions:

```
import utils

data = get_data()
metadata = {
    'capname': 'Example Capture',
    'lasernum': 3,
    'lasername': 'Clara'
}

utils.save("outfile", data, metadata)

# Then later, when you need the data back:
```

```
data, md = utils.load("outfile")
```

3.2 Editing Files

The script `edit_md.py` can be used to modify a file's settings. Run it as

```
python scripts/edit_md.py filename
```

and follow the instructions on the screen.

4 Settings

Various settings for the project can be found in `settings.py`. The first three entries in the file (`ROOT_DIR`, `FIGS_DIR`, and `DATA_DIR`) set the working directory, and directories to save and load figures and data from, respectively. The snippet of Python code for `ROOT_DIR`

```
ROOT_DIR = Path(__file__).parents[1]
```

simply sets the root directory to the parent of the scripts directory. You shouldn't need to change any of these parameters, short of a major reorganization of the project.

The rest of the settings are snippets from the VISA ids of various instruments used:

- The scope is a Tektronix TDS-2000C 4-channel scope w/ USB out. There's one (and only one) of these in the JLab.
- The function generator is an Agilent 33210A.
- The multimeter is an

I believe these snippets will be the same for any instrument of the same type, but I am not sure. If you get "InstrumentNotFoundError" or similar, this may be the cause.

5 Characterizing MZMs

To characterize an MZM, connect it to the associated laser, and turn on the laser power. Feed a $10V_{pp}$ triangle wave into the MZM at 8Hz. Connect channel 2 on the scope to the function generator, and channel 1 to the detector output. Run

```
python scripts/MZ_getdata.py
```

and follow the prompts for filename and power level.

After collecting this capture, run

```
python scripts/MZ_characterize.py <filename>
```

to fit the appropriate curve. The fitted parameters will be printed to the terminal, and the plot saved to the figures directory. Make sure to inspect the plot to verify that the fit actually worked!

The number you probably most care about is the $V\pi$ figure reported on the terminal—this can be compared to the value on the datasheet, and used to determine the necessary offset voltages.

6 Characterizing Lasers

Unfortunately, this is the least-automated procedure in the whole project. Because I wasn't able to find out how to programmatically set the laser controller (I believe it's possible, but would require an Arduino or something), it requires manually measuring voltage at a variety of input currents. That being said, the procedure is:

1. Connect the output of the laser to the detector, without the polarization controller or MZM in the fiber path.
2. Connect the output of the detector to a multimeter set to measure DC voltage.
3. Open a spreadsheet, and create two columns (current and voltage).
4. At various input currents, measure the output voltage
 - (a) There is a threshold current, typically about 8mA
 - (b) My measurements were spaced at 1mA, but decent results should be achievable with slightly wider spacing
 - (c) There is an overcurrent protection around 50mA, but you shouldn't need to go that high anyways—typically, current while running experiments is ca. 20mA.
5. Save the spreadsheet (as a CSV) in `data/lasers`
6. Create a `.dat` file for the laser and detector configuration. Entries in this file:
 - (a) `name`: The nickname given to the laser (Alice, Clara, and Kira)
 - (b) `wavelength`: The peak wavelength, as given by the spec sheet
 - (c) `r`: The detector efficiency at that wavelength, as given by the detector spec sheet
 - (d) `g`: The gain setting of the detector
 - (e) `threshold`: The expected threshold current, from the datasheet.
 - (f) `opCurrent`: The expected current-power relation, in mW/40mA. Also from the datasheet.

Use one of the existing files as a formatting guide.

7 Characterizing Symetrix 402 Delay Module

For the delay module, we want to know whether and how the input signal is distorted and shifted (apart from the time delay, obviously). Some tests, such as determining the maximum input voltage and observing distortion of non-sinoidal signals, are purely qualitative, and need little explanation.

Quantitative tests we can perform on the delay module include measuring its gain and phase shift across a variety of frequencies. To run these tests, connect a function generator to the input of the delay module, and connect the output (I used output 1, but it should not matter) to a scope. Also run the signal directly from the function generator to the scope. Connect both the function generator and scope to the computer via USB, and run

```
python scripts/freqsweep.py <fmin> <fmax> <df> <filename>
```

replacing each of the arguments as appropriate. What this does is, for each frequency in

$$[f_{\min}, f_{\min} + df, f_{\min} + 2df, \dots, f_{\max}],$$

generates a sine wave of that frequency and saves the resulting waveforms. Because the entire waveform is saved for each frequency, the data from these captures tend to be somewhat large.

Once the data have been collected, we can plot the relationships between frequency and gain and phase. The script `delay_freq_phase_diff.py` determines the phase shift induced in the signal as a function of frequency, and the script `delay_gain.py` determines the overall gain of the circuit, as well as the probability that the gain is frequency-independent (tests show that, although the effect of frequency is minuscule, it does slightly affect the gain).

8 Characterizing bandpass filter

To determine the parameters of the bandpass filter, we need to create a Bode plot, plotting gain against frequency. Run

```
python scripts/freqsweep.py
```

as described in section 7 above. Then run

```
python scripts/filter_gain.py
```

on the resulting file. With the bandpass filters in this circuit, centered around 8Hz, a frequency range of 4-12 Hz should be sufficient.

9 The Circuit

Figure 2 shows the full form of the circuit. The left and right sides are identical, each being the circuit for one of the loops. The center contains the AD5206 digital potentiometer IC, as well as the Arduino used to drive it. The (many, many) intricacies of the AD5206 are detailed in a later section. Other than that, though, the circuit is relatively straightforward. All 8-DIP chips are standard 411 op-amps, and other than that, all components are passive. Going top to bottom, the components of each loop are:

- Op-amp 1: Voltage follower for digipot.
- Op-amp 2: Voltage follower for digipot
- Op-amp 3: Summing amplifier, with maximum gain of 4. The resistors from the voltage followers are both $27\text{ k}\Omega$, while the self-feedback resistor is $100\text{ k}\Omega$.
- Op-amp 4: Bandpass filter, with nominal parameters $f_0 = 8\text{ Hz}$, $A_0 = 1.0$, $Q = 4.5$. The filter design is the multiple-feedback active bandpass filter described in section 11.6.1 of [1], and is shown in figure 3. Note that, in the actual circuit, many of the resistors are actually two resistors in series, to allow for a closer approximation of the desired resistance.
- Op-amp 5: Unity-gain summing amplifier, for introducing a DC offset. All resistors are $82\text{ k}\Omega$.

9.1 The AD5206 Digital Potentiometer


The AD5206 digital potentiometer is—by far—the most finicky part of the circuit. I accidentally fried a number of these chips while trying to determine how they worked, and I know Yunjia faced similar problems. A brief summary application note on these devices I wrote earlier in the summer:

- The AD5206 behaves like six mechanical potentiometers.
- Each potentiometer has three terminals (A,B, and W), just like a mechanical.
- These potentiometers are (to the best of my knowledge) completely isolated—I did not observe any leakage between them.
- There are two power terminals on the chip: V_{SS} and V_{DD} .
- These terminals are subject to the constraints that

$$V_{SS} \leq 0 \leq V_{DD}$$


and

$$|V_{DD}| + |V_{SS}| \leq 5.5V$$



fritzing

Figure 2: The electronic part of the optoelectronic oscillator



fritzing

Figure 3: The bandpass filter, with (nominal) component values

- The potentials of all three wiper terminals MUST be between V_{SS} and V_{DD} :

$$V_{SS} \leq \{A, B, W\} \leq V_{DD}$$

- If this condition is not met, the response is nonlinear, and signals are distorted.

Monopolar Mode

These potentiometers are usually used in (and most online resources assume) monopolar mode

- In this mode, $V_{SS} = 0V = \text{GND}$, and $V_{DD} = 5V$.
- This is very convenient, since the logic output of most microcontrollers (e.g. Arduino) is 5V, so no step-down is needed.
- However, because of the voltage bounds mentioned above, this means we can't send zero-offset AC signals (like we're working with here) through the digipot
- We could add an offset using an opamp, but that makes scaling very difficult, and adds unnecessary noise.

Bipolar Mode

Instead, we can operate in bipolar mode, where $V_{SS} < 0$. Some considerations in this mode:

- The potential bounds still apply. If we power the digipot with $\pm 2.5V$, the signal must fall within those bounds as well. This should not be a problem, since our signal should be within the $\pm 1V$ envelope.
- The logic voltage must be within $0.3V$ of V_{DD} (or specifically, when $V_{DD} = 3V$, the logic must be between $2.6V$ and $3.3V$, as per the datasheet).
- Logic low is still ground.
- I am driving the chip with an Adafruit ItstBitsy arduino, which uses 3.3V logic. I have implemented voltage dividers on each data channel to bring the signal to the required 2.5V.

Communication

- We communicate with the potentiometer using SPI (serial peripheral interface).
- To set a resistance:
 - Drive the CS (chip select) pin low
 - Write the address as a single byte, MSB first. Because valid addresses are 0 – 5 (zero-indexed), the word will look like `0b00000101` (for address 5).
 - In Arduino, this can be done with `SPI.transfer(channel)`, where channel is in the range 0 – 5.
 - Then write the desired potentiometer value as a single byte, ranging from 0 (minimum resistance) to 255 (maximum resistance). In theory, the output resistance at step n will be

$$R_{\text{out}} = (n/255) \cdot 50\text{k}\Omega$$

Experimental response curves are detailed below.

- Drive the CS pin high.
- The code used on the Arduino is in the associated Github repository.

Characterizing and Using the Chip

- Be careful when testing the chip. A natural instinct would be to test the chip by hooking it up to a multimeter and sweeping through the various settings. However, multimeters (typically) work by applying a 5V potential difference, and measuring the associated current. Since (in bipolar mode) 5V is outside of the $[V_{ss}, V_{dd}]$ range, this will kill the chip. I fried at least two chips doing this.
- There are two ways we can use the chip. In rheostat mode, we use the chip directly as a variable resistor that can be digitally tuned, such as by using it as one of the resistors in a summing amplifier. In the potentiometer mode, we use the chip as a programmable voltage divider—essentially, a variable gain amplifier with a gain between 0 and 1.
 - In rheostat mode, we need to worry about overcurrent. Assuming a 1V signal (broadly typical), and that the output is tied to ground (which, with the summing amplifier design, is true), the current flowing through the chip at the minimum resistance of $50\,\Omega$ is

$$I = \frac{V}{R} = 0.02\,\text{A} = 20\,\text{mA}$$

Given that the absolute maximum continuous current should not exceed 2.5 mA, this is bad. I killed a chip this way too (it got very warm).

- In potentiometer mode, terminal B is tied to ground, and terminal A is tied to V_{in} . The output voltage on the wiper is then

$$V_{out} = \frac{n}{255} V_{in}$$

where n is the digital input.

The other main advantage of potentiometer mode is reduced temperature dependence. The digital potentiometer (and all potentiometers) has a relatively large temperature coefficient, so the actual resistance/tick value changes with temperature. However, in the voltage divider equation, multiplying all resistances by some factor does not change the output voltage—so, by using this mode, we do not have to account for temperature in our circuit.

Response curves

Potentiometer	Base resistance (Ω)	Response (Ω/tick)
Ideal	50	196.1
R1	30.3	202.6
R2	39.5	202.3
R3	39.9	202.1
R4	42.0	202.1
R5	10.9	202.6
R6	27.4	202.4

These resistors have a small (near-negligible) wiper resistance, and a slightly higher maximum resistance than advertised. The typical maximum resistance was about 51.5 k Ω . No statistically significant nonlinearity was observed.

Other notes

- It may be worthwhile to look into other digital potentiometers, especially if this one breaks or burns out. As far as I can tell, Analog Devices (AD) is the largest and most reputable manufacturer of these chips.

One chip I was looking at was the AD5263, which only has four channels (instead of the six found on the AD5206), but can be powered by a $\pm 5V$ supply (and therefore 5V logic).

- Unfortunately, Analog seem to have discontinued all DIP-packaged chips, so any new model of potentiometer (or replacements for the existing one, if it breaks) means dealing with SMD parts and breakout boards.
- It may be worthwhile to look into replacing the digital potentiometer with variable gain amplifiers, especially since the potentiometer/voltage follower combo essentially act as a variable gain amplifier anyways.

10 Utilities

10.1 `combine_figs.py`

This is a useful little script for combining multiple figures together on one page for printing or easy comparison. To use:

```
python scripts/combine_figs.py <outfile.pdf> <fig1> <fig2> ...
```

where all the figures are in the same directory as `outfile.pdf`. This thing uses LaTeX behind the scenes, so if things go south, that's probably the place to look.

References

- [1] Robert E. Simpson. *Introductory Electronics for Scientists and Engineers*. Addison-Wesley, 1987.